



GE Medical Systems
Information Technologies

CARDIOSOFT THE ECG FOR YOUR PC

Feature Highlights

- Easy installation
- Designed for the future
- Can be integrated into office programs and hospital information systems
- Fully WEB-enabled



All you need for a complete ECG workstation is a PC and CardioSoft. CardioSoft consists of the application software and the CORINA ECG amplifier. The installation is straightforward and easy because the ECG amplifier connects to a standard parallel port. There is no need to make changes to the PC hardware.


A variety of software options enable the system to be configured individually for any requirements. Because of the graphical user interface under Microsoft Windows, you will find the system readily accessible and easy to use.

All data is presented quickly and clearly on one screen. Events are indicated in colour.

DDE and OLE functions facilitate easy data exchange with other Windows programs.

The pump for the KISS electrode application system can optionally be placed into the ECG amplifier's case in order to save space and time.

Highlights



CARDIOSOFT THE ECG FOR YOUR PC

CardioSoft is tailored to your needs and can be integrated into existing and future systems.

The ECG system can be connected to a network and is modular, so it can grow with your requirements. You can have access to data and values at any time from any place. The management, documentation, revision, evaluation and annotation of results is simplified immensely, and providing more effective and accurate results.



Furthermore, the integrated GDT/BDT interface offers you the possibility of exchanging data with practice programs and hospital information systems. Multiple data entry can be eliminated and the amount of work reduced considerably.

Transfer the most important data for further evaluation into word processing or spreadsheet programmes using the export function. Thus an economic and effective management of all administrative processes in the clinic and practice is achieved.

Integration into a network is very easy. The system is compatible with standard hardware and software structures.

The WEB option gives you access to the complete stored data by using any standard internet browser. With this platform independent technology you can display reports as well as waveforms from any place at any time.

By using CardioSoft today your work will not only be considerably simpler and safer, but will also offer the best prospect for effective use in the future.



GE Medical Systems
Information Technologies

gemedicalsystems.com

European Headquarters
GE Medical Systems
Information Technologies GmbH
P.O. Box 60 02 65
79032 Freiburg • Germany
Tel. +49 761 45 43 - 0
Fax +49 761 45 43 - 233

World Headquarters
GE Medical Systems
Information Technologies, Inc.
8200 West Tower Avenue
Milwaukee, WI 53223 • USA
Tel. +1 414 355 5000
Fax +1 414 355 3790

Asia Pacific
GE Medical Systems Hong Kong Ltd.
11th Floor, The Lee Gardens, 33 Hysan Ave.
Causeway Bay • Hong Kong
Tel. +852 2100 6300
Fax +852 2100 6292